

# Swabs: types and codes

**Patient Request Forms** and **Swabs** should be labelled with the body site from which the sample was taken. **This is important**. The swab site determines the appropriate culture media required to target the most likely pathogens.

SITE CODE SAMPLETYPE

#### **Culture Swabs**

Candida only swab	CANC	Black or Blue Micro Swab
Cervical swab	CERS	Black or Blue Micro Swab
Ear swab	EARS	Blue or Orange Micro Swab
Eye swab	EYES	Blue or Orange Micro Swab
Gonorrhoea	GONN	Black Charcoal Swab
High vaginal swab	HVS	Black or Blue Micro Swab
Nasal swab	NASS	Blue or Orange Micro Swab
Oral swab	ORSW	Black or Blue Micro Swab
Penile swab	PENS	Black or Orange Micro Swab
Rectal swab	RECG	Black or Blue Micro Swab
Skin swab	SKIS	Blue Micro Swab
Throat swab	THRS	Blue Micro Swab
Urethral swab	URES	Black or Orange Micro Swab
Vaginal swab	VAGS	Black or Blue Micro Swab
Vulval swab	VULV	Black or Blue Micro Swab
Wound swab	WOUS	Black or Blue Micro Swab

### **MRSA by Culture**

MRSW	<b>Blue</b> Micro Swab x 1 - state site
MRW2	<b>Blue</b> Micro Swab x 2 - state sites
MRW3	<b>Blue</b> Micro Swab x 3 - state sites
MRW4	Blue Micro Swab x 4 - state sites
MRW5	<b>Blue</b> Micro Swab x 5 - state sites

### **Rapid MRSA by PCR**

MRSA	<b>Blue</b> Micro Swab x 1 - state site
MRS2	<b>Blue</b> Micro Swab x 2 - state sites
MRS3	<b>Blue</b> Micro Swab x 3 - state sites
MRS4	Blue Micro Swab x 4 - state sites
MRS5	<b>Blue</b> Micro Swab x 5 - state sites

#### **Blue Micro/Transwab**

are multipurpose, culture swabs in transport medium

# Orange Micro/Transwab

are small, thin wire culture swabs in transport medium

**Black Charcoal Micro/Transwab** Wound, skin and urogenital.